CE₀₄₂₆

Aghi dentali universali

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CE₀₄₂₆

Dental universal needles

Needles for dental use with hub treading, european standard for tube-ampule syringe, sterilized by gamma irradiation, silicone treated, disposable. The welding point which joins the safety cap to the protector guarantees that the needle remains sterile.

INDICATIONS

These needles must only be used by dental surgeons for the administration of local dental anaesthesia.

INSTRUCTION FOR USE

Positioning of the needle on the syringe:

- 1) Break the seal
- 2) Twist and pull to remove the safety cap
- 3) Insert the needle into the syringe, being careful not to remove the protector.
- 4) Remove the protector by twisting gently, just before the administration of the anaesthesia.

AFTER USE:

- 1) Take the needle out of the syringe very carefully, in order to avoid pricking anyone accidentally.
- 2) Do not leave the needle lying around, but dispose of it immediately in a suitable container which conforms to regulations

CAUTION

Do not use if the protection seal is broken. Needles after use carry a high risk of contamination; be extremely careful not to prick anyone accidentally. These needles must be used once only and cannot be resterilized.

WARNINGS

The safety cap must never be replaced on the needle by hand, or there could be the risk of accidentally pricking someone. Once used, the safety cap must not be replaced on the needle, but the needle must be disposed of in a suitable container which conforms to regulations, without bending, breaking or touching it in any other way. Conserve at room temperature, away from sources of direct heat and damp places. The needle is used for the administration of dental anaesthetic.

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