

MATERIAL SAFETY DATA SHEET

Tubulitec Cavity Liner

Date of issue: 2004-09-13

Page1/3
Replaces date: 2003-05-08

1. Identification of the substance/preparation and of the company/undertaking

PRODUCT NAME:	Tubulitec Cavity liner
TYPE OF PRODUCT:	Dental Product. For insulation of dental cavities.
MANUFACTURER:	Dental Therapeutics AB
ADDRESS:	Paviljongvägen 3-5
POSTAL ADDRESS:	SE-132 40 Saltsjö-Boo
TELEPHONE:	+46 8 7168141
TELEFAX:	+46 8 7182175
INFORMATION ISSUED BY/CONTACT PERSON:	Kenneth Larsson

2. Composition /Information on ingredients

No:	Name:	CAS-nr:	Conc:	Classification:	R-phrases
1	Ethyl acetate (C ₄ H ₈ O ₂)	141-78-6	80%	F	11
2	Calcium Hydroxide	1305-62-0	5%	Xi	41
4	Calcium oxide	1305-78-8	1%	Xi	41
5	Diiiodithymol	552-22-7	1%	Xi	36/38

T+= Very toxic. T=Toxic. C= Corrosive, Xn= Harmful, Xi= Irritating, E= Explosive, O= Oxidising, F+= Extremely flammable, F= Highly Flammable.

3. Hazards Identification

HIGHLY FLAMMABLE.

RISK OF SERIOUS EYE DAMAGE.

IRRITATING TO SKIN.

4. First aid measures

After inhalation: Fresh air. If symptoms remains seek medical advice.

After skin contact: Wash of with soap and water.

After ingestion: Avoid vomiting. Make victim drink lemon juice, table vinegar, or raw eggs mixed into milk, followed by plenty of water. Seek medical advice.

After eye contact: Rinse out with plenty of water for at least 10 minutes with the eyelid held wide open. If symptoms remains, seek medical advice.

5. Fire-fighting measures

Suitable extinguishing media: Powder, foam.

Extinguishant: Carbon dioxide, powder for liquid.

Special risk: Combustible. Vapour heavier than air. Formation of explosive mixtures possible with air. Keep away from sources of ignition. Take measures to prevent electrostatic charging.

6. Accidental release measures

Person-related precautionary measures: Do not inhale vapours. Room temperature vapours can form explosive mixtures with air. Be aware of the risk of explosion.

Procedures for cleaning/absorption: Take up with liquid absorbent material. Forward for disposal. Clean up.

7. Handling and storage

Handling: No further requirements.

Storage: Tightly closed. In well-ventilated place. Keep away from sources of ignition and heat. At room temperature.

MATERIAL SAFETY DATA SHEET

Tubulitec Cavity Liner

Date of issue: 2004-09-13

Page 2/3
Replaces date: 2003-05-08

8. Exposure controls/personal protection

Name:	CAS-No:	Max allowable conc. :
Ethyl acetate	141-78-6	1400 mg/m ³
Zinc Oxide	1314-13-2	5mg/m ³
Calcium Oxide	1305-78-8	5mg /m ³

Specific control parameter: Possibilities to eye flush near the work place.

Eye protection: Required

Hand Protection: Required

Skin Protection: Wear suitable protecting clothes.

9. Physical and chemical properties

Form: Liquid	
Colour: Yellow-brown	
Odour: Pleasant fruit-like	
Melting temperature: -83 degr. C°	
Boiling temperature: +77 degr. C°	
Density:	
Flash point: -5 degr.C°	
Explosion limits: %--%: 2,1-11,5 vol%	
pH value: 7,3	
Vapor pressure: 9,7 kPa	
Solubility in Water:	

10. Stability and reactivity

Conditions to be avoided: None

Substances to be avoided: Alkali metals, hydrides, strong oxidizing agents, water with air and light.

Hazardous decomposition products: Ethanol, acetic acid

11. Toxicological information

Acute toxicity:

Ethyl Acetate:	LD50 (oral,rat): 5620 mg/kg
Calcium Hydroxide:	LD50 (oral, rat): 7340 mg/kg
Diiodidtymol:	LD50 (oral, rat): 950 g/kg

Further toxicological information:

After ingestion and inhalation: Mucosal irritations, lack of appetite, headache. In high concentrations: salivation, nausea, vomiting, narcosis, respiratory paralysis.

After skin contact: Degreasing effects on the skin, possibly followed by secondary inflammation. After long-term exposure: sensitization with allergic manifestations.

12. Ecological information

Biologic degradation:

Biologic degradability: Good

Exotoxic effects: When used properly, no impairments in the function of waste-water-treatment plants are to be expected.

Toxic for aquatic organisms: Fish: LC50:333 mg/l.

Further ecologic data: Product reacts with water: fish EC0:50 mg/l (acetic acid)

Do not allow entering drinking water supplies, waste water, or soil!

13. Disposal considerations

Contact waste disposal companies that will advise you how to dispose of special waste.

MATERIAL SAFETY DATA SHEET

Tubulitec Cavity Liner

Date of issue: 2004-09-13

Page 3/3
Replaces date: 2003-05-08

14. Transport information

ADR/RID (land)
UN no: 1173
Class: 3, 3b
AND/ADNR (River)
Not examined
IMGD (Sea)
UN no: 1173
Class: 3.2
Packaging group: II
ICAO-IIATA and IATA-DGR (air)
UN no: 1173
Class: 3
Packaging group: II

15. Regulatory information

Labelling according to Council Directive 93/42/EEC Annex V and the standards EN ISO 9002 and EN 46002.

Symbol: F, Xi

Classification: Highly flammable, Irritating.

R-phrases: 11, 41, 36/38 Highly flammable, Risk of serious damage to eye. Irritating to eyes and skin.

S-phrases: 16-23-29-33 Keep away from sources of ignition – No smoking. Do not breathe vapor. Do not empty into drains. Take precautionary measures against static discharges.

Ethyl Acetate EEC-No: 205-500-4
CAS-No: 141-78-6

THE BOTTLE CONTAINS 10 ML

16. Other information

ONLY FOR PROFESSIONAL DENTAL USE

The glass container and the cap is both recyclable.

The information above is believed to be accurate and represents the best information currently available to us. Users should make their own investigations to determine the suitability of the information for their particular purposes. In no way shall the company be liable for any claims, losses, or damages of any third party or for lost profits or any special, indirect, incidental, consequential or exemplary damages.